



蘇州恒祥進出口有限公司

SUZHOU HENGXIANG IMPORT & EXPORT CO., LTD.

**ONE STOP MEDICAL PRODUCTS SUPPLY**

- Aiming at supplying one-stop professional service to you for medical products, especially for medical disposables like:

*Disposable Medical Device...*  
*Disposable Medical Dressings...*  
*Disposable Medical Non-woven Products...*  
*Disposable Medical Gloves...*

- Afford different payment terms, especially OA term (Open Account Term) for supporting your cash flow by maximum, as long as you have good credit. The longest time could be OA 90days
- Specially good at handling different medical products from different factories but ship in one container
- High quality but cheap price medical products with fast delivery time  
All selected factories (our suppliers) are around our cities (in 1~3 hours by car, these factories are equipped with modern production machines under modern management. All these factories have secured CE, SO, China FDA, USA FDA etc.)
- All ordered medical products would be strictly supervised by us, whether in the course of "before production". "ongoing production or "finished products"
- Over 40 years import and export business in medical field (used to be a state-run company which specialized in import & export of medical product 40 years ago, ended till year 2001)
- A member of China Chamber of Commerce for Import & Export of Medicines & Health Products
- Turnover of year 2023 reaches 150 million USD
- Locates famous Suzhou city, only 80km far from Shanghai (30 minutes by High-speed train or 1 hour by car) Suzhou is an old but most beautiful city with 2500 years history, used to be ancient China's economy center, culture center and art center!  
But now Suzhou city is the China No. 1 or 2 industrial city too.
- Registered medical brand "HEPPO" is already a famous in some countries and regions



Wechat:  
008613063784080



Whatsapp:  
008613063784080



**HEPPO**<sup>®</sup>  
ONE STOP MEDICAL PRODUCTS SUPPLY

蘇州恒祥進出口有限公司  
SUZHOU HENGXIANG IMPORT & EXPORT CO., LTD.

Tel: 0086-512-68626864  
Fax: 0086-512-68621176  
Mob: 0086-130-63784080  
E-mail: pzw@hengxiang.cn  
panjacky1979@163.com  
Website: www.heppo.cn



潘振偉 | Jacky Pan  
銷售經理 | Sales Manager

(wechat) +8613063784080  
(QQ) 81744697  
(skype) jackypan1979  
(whatsapp) +8613063784080

江蘇省蘇州市吳中區澄湖路888號A座10樓, 11樓  
Add:10-11F, Building A, No.888 Chenghu Road, Wuzhong District, Suzhou, Jiangsu, China

郵編/Zip Code: 215168

for more catalogues & price list, please scan   
or visit: <http://www.hengxiang.cn/en/file-80793.html>




For Catalogues & Prices

## Disposable Hypodermic Disposable Syringe Needle

### Sizes & Specifications:


Size		Color	Regular Needle Length						
OD	Gauge		1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"	
(mm)			13mm	16mm	19mm	25mm	32mm	38mm	
0.30	30G	Light Yellow	•						
0.33	29G	Red	•						
0.36	28G	Blue-Green	•						
0.40	27G	Grey	•						
0.45	26G	Brown		•					
0.50	25G	Orange		•	•				
0.55	24G	Purple			•	•			
0.60	23G	Blue				•	•		
0.70	22G	Black					•	•	
0.80	21G	Green						•	•
0.90	20G	Yellow						•	•
1.10	19G	Beige						•	•
1.20	18G	Pink						•	•
1.60	16G	White						•	•

- Packing: 1pc/soft blister bag, 100pcs/box, 5000pcs/ctn
- Sterile: EO Sterile
- Needle for animal use is also available




## Disposable 3-Part Syringe

- Sizes: 1ml, 2ml, 2.5ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml, 60ml
- Syringe Tip: Luer Tip / Luer Lock / Catheter Tip
- Needle: with / without Needle
- Piston: Latex / Latex-Free
- Plunger: White / Semi-transparent
- Individual Package: 1pc/PE bag or soft blister bag
- Mid-Packing: Papercard Box or Big PE Bag
- Sterile: EO Sterile
- Remarks:
  - a) For 5ml~60ml, eccentric tip is available
  - b) For 1ml Syringe, tuberculin syringe and BCG vaccine are also available



## Disposable 2-Part Syringe


- Latex Free
- Sizes: 2ml, 2.5ml (Centric Tip) ; 5ml, 10ml, 20ml (Eccentric Tip)
- Syringe Tip: Luer Tip
- All others are same as Disposable 3-Part Syringes



A: 2ml/2.5ml (Centric Tip)  
B: 5ml (Eccentric Tip)  
C: 10ml (Eccentric Tip)  
D: 20ml (Eccentric Tip)



## Disposable Insulin Syringe

- Ultra fine needle for minimizing patient's pain
- Accurate injection qty
- With normal detachable needle (**Model A**) or with fixed & jointed needle (**Model B**)
- Sizes: 0.5ml, 1ml with needle 29Gx1/2" or 30Gx1/2" (40unit, 80unit, 100unit)
- All others are same as Disposable 3-Part Syringe



## Disposable Syringe & Needle for Medicine-Solving

- For Hand Use: 5ml, 10ml, 20ml, 30ml, 50ml, 100ml with Special Needle 16G or 18G
- For Machine Use: 20ml with Special Needle 16G or 18G



Needle for Medicine-Solving



Needle for Animal Use



# Infusion Set with Burette





**B-100A**



**B-100B**



**B-150A**



**B-150B**



**B-150C**

## Configuration of Burette Infusion Set

**Remarks:**

**1) Standard Configuration (Economic Configuration):**

spike, hanging ropes, drip chamber (50ml /100ml/150ml, 60drops/ml), tubing (length: 25+150+30cm, O.D/I.D: 4mm/3mm), two regulators (50mm), male luer-lock, with DEHP ingredient

**2) Optional Components to Add:** T-shape spike, air inlet, regulator (50mm/55mm with tube clasp); floating disk valve (latexfree isoprene rubber /latexfree silicone rubber), solution filtration disk; Y-shape injection site(latexfree); needlefree injection site(latexfree); extension tubing (O.D/I.D: 4mm/3mm); frosted elastic tubing; with or without "hypodermic needle, scalp vein set, anti-reflux valve, blue vented cap, micro flow regulator, micro filter" etc; light shielding function etc.; DEHP-free

Model No.	Code of Configuration	Nutral Latex or Latexfree	Chamber Capacity	Spike	Hanging Rope	Upper Regulator (Roller Clamp)	Floating Disk Valve	Drip Chamber	Solution Filter Disk	Under Regulator (Roller Clamp)	Y-shape Injection Site	Adaptor	Needle	Remarks
B-50-1F	JR1-50F2J1YTOV	latexfree	50ml	w/o air inlet	with	50mm	silicone rubber	standard 60drops/ml	/	55mm with tube clasp	with	luer tip	w/o needle	with scalp vein set
B100-1F (B-100A)	AR1-100F1J1YLO	latexfree	100ml	with air inlet	with	50mm	isoprene rubber	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	
B150-2F	AR1-150F0J1LOC2	latexfree	150ml	with air inlet	with	50mm	/	standard 60drops/ml	with	55mm with tube clasp	/	luer lock	w/o needle	with 2 needlefree injection sites
B150-3F	AR1-150F2YLO	latexfree	150ml	with air inlet	with	50mm	silicone rubber	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	
B150-4F	AR1-150F0J1YLO	latexfree	150ml	with air inlet	with	50mm	/	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	chamber with blue-strip
B150-5F	AR1-150F0J1YLO	latexfree	150ml	with air inlet	with	50mm	/	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	for blood transfusion!
B150-6F	AR1-150F0J1YLOZ	latexfree	150ml	with air inlet	with	50mm	/	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	light-shielding
B150-7F	AR1-150F0TO	latexfree	150ml	with air inlet	with	50mm	/	/	/	/	/	tube	w/o needle	
B-100A (B100-1F)	AR1-100F1J1YLO	latexfree	100ml	with air inlet	with	50mm	isoprene rubber	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	
B-100B	AR1-100F2J1YLO	latexfree	100ml	with air inlet	with	50mm	silicone rubber	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	
B-150A	AR1-150F0J1YLO	latexfree	150ml	with air inlet	with	50mm	/	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	
B-150B	AB-150F2J1YLO	latexfree	150ml	with air inlet	with	Robert Clamp	silicone rubber	standard 60drops/ml	with	55mm with tube clasp	with	luer lock	w/o needle	
B-150C	AJR1-150F0TO	latexfree	150ml	T-spike with air inlet	with	50mm	/	/	/	/	/	tube	w/o needle	



# Infusion Set





54-1F



54-2F



54-3F



54-4F



54-5F



54-6F





54-7F



54-8F



55-3F



62-1F



62-2F



62-3F









## Infusion Set




## Enteral Feeding Set













# Componnts of Spike, Drip Chamber and Regulator

## Spike & Drip Chambers



- A-drip chamber 50mm (with air inlet)
- B-drip chamber 50mm (w/o air inlet)
- C-drip chamber 50mm (with air inlet)
- D-drip chamber 54mm (with air inlet, with wings)
- E-drip chamber 54mm (w/o air inlet, with wings)
- F-drip chamber 72mm (with air inlet, with wings)
- G-drip chamber 62mm (with air inlet)
- H-drip chamber 72mm (with air inlet)
- J-drip chamber 72mm (w/o air inlet)
- K-drip chamber 72mm (with air inlet, light shielding)
- L-drip chamber 89mm (with air inlet, with ring)

## Regulators (Roller Clamps)



- R1-50mm
- J1-55mm with tube clasp
- R2-55mm safety



## Configuration of Infusion Set

### Remarks:

1) **Standard Configuration (Economic Configuration):** spike, drip chamber (50mm, 20drops/ml), tubing (length: 150cm, O.D/I.D: 4mm/3mm), regulator (50mm), luer-tip adaptor, with DEHP ingredient

2) **Optional Components to Add:** T-shape spike, spike with wings; air inlet; drip chamber (54mm/62mm/72mm/89mm. 10drops/ml, 60drop/ml); regulator (50mm/55mm with tube clasp/ 55mm with safety hole); solution filtration disk;natural latex /isoprene rubber (latexfree) injection bulb (line-shape/pear-shape); Y-shape injection site(latexfree); needlefree injection site(latexfree); 3-way stopcock injection site(latexfree); extension tubing (O.D/I.D: 4mm/3mm); frosted elastic tubing; male luer lock adaptor; with or without "hypodermic needle, scalp vein set, anti-reflux valve, blue vented cap, micro flow regulator, micro filter" etc; light shielding function etc.; DEHP-free

Model No.	Code of Configuration	Nutral Latex or Latexfree	Spike	Drip Chamber	Solution Filtration Disk	Regulator (Roller Clamp)	Natural Latex Injection Site	Y-shape Injection Site	Adaptor	Needle	Other Components	Remarks
50-1F (2#)	A50GR1PTW	latex	with air inlet	50mm	with	50mm	pear-shape		luer tip	with needle		
50-2F (6#)	A50GR1YWLW	latexfree	with air inlet	50mm	with	50mm		with	luer lock	with needle		
50-3F (3#)	A50GR1PLW	latex	with air inlet	50mm	with	50mm	pear-shape		luer lock	with needle		
50-4F	50GR1YWLW	latexfree	w/o air inlet	50mm	with	50mm		with	luer lock	with needle		
50-R	50J1PTW	latex	w/o air inlet	50mm		55mm with tube clasp	pear-shape		luer tip	with needle		Indonesia/Japan Model
50-H	W54J1STO	latexfree	w/o air inlet with wings	54mm		55mm with tube clasp	pear-shape		luer tip	w/o needle		isoprene rubber injection bulb
54-1F	AW54GR1PTW	latex	with air inlet with wings	54mm	with	50mm	pear-shape		luer tip	with needle		
54-2F (12#)	AW54GJ1YLO	latexfree	with air inlet with wings	54mm	with	55mm with tube clasp		with	luer lock	w/o needle		
54-3F	AW54J1SYTO	latexfree	with air inlet with wings	54mm		55mm with tube clasp	pear-shape	with	luer tip	w/o needle		Y & isoprene rubber injection site
54-4F	54GJ1YLO	latexfree	w/o air inlet	54mm	with	55mm with tube clasp		with	luer lock	w/o needle		
54-5F	AW54GJ1YLOF	latexfree	with air inlet with wings	54mm	with	55mm with tube clasp		with	luer lock	w/o needle	with mirco flow regulator	
54-6F	AW54J1STOM	latexfree	with air inlet with wings	54mm		55mm with tube clasp	pear-shape		luer tip	w/o needle	with micro filter	isoprene rubber injection bulb
54-7F	AW54GJ1LO3	latexfree	with air inlet with wings	54mm	with	55mm with tube clasp			luer lock	w/o needle	with 3-way stopcock	
54-8F	AW54GJ1YLOZ	latexfree	with air inlet with wings	54mm	with	55mm with tube clasp		with	luer lock	w/o needle		light shielding
55-1F (F)	A89GJ1YLO	latexfree	with air inlet	89mm (with ring)	with	55mm with tube clasp		with	luer lock	w/o needle		two-part chamber with rings
62-1F	AW62GJ1LO	latexfree	with air inlet with wings	62mm	with	55mm with tube clasp			luer lock	w/o needle		no injection site

## Configuration of Infusion Set

### Remarks:

1) **Standard Configuration (Economic Configuration):** spike, drip chamber (50mm, 20drops/ml), tubing (length: 150cm, O.D/I.D: 4mm/3mm), regulator (50mm), luer-tip adaptor, with DEHP ingredient

2) **Optional Components to Add:** T-shape spike, spike with wings; air inlet; drip chamber (54mm/62mm/72mm/89mm. 10drops/ml, 60drop/ml); regulator (50mm/55mm with tube clasp/ 55mm with safety hole); solution filtration disk;natural latex /isoprene rubber (latexfree) injection bulb (line-shape/pear-shape); Y-shape injection site(latexfree); needlefree injection site(latexfree); 3-way stopcock injection site(latexfree); extension tubing (O.D/I.D: 4mm/3mm); frosted elastic tubing; male luer lock adaptor; with or without "hypodermic needle, scalp vein set, anti-reflux valve, blue vented cap, micro flow regulator, micro filter" etc; light shielding function etc.; DEHP-free

Model No.	Code of Configuration	Nutral Latex or Latexfree	Spike	Drip Chamber	Solution Filtration Disk	Regulator (Roller Clamp)	Natural Latex Injection Site	Y-shape Injection Site	Adaptor	Needle	Other Components	Remarks
62-2F	AW62GJ1SLO	latexfree	with air inlet with wings	62mm	with	55mm with tube clasp	pear-shape		luer lock	w/o needle		
62-3F	AW62GJ1SYLO	latexfree	with air inlet with wings	62mm	with	55mm with tube clasp	pear-shape	with	luer lock	w/o needle		Y & isoprene rubber injection site
62-4F	AW62GJ1Q	latexfree	with air inlet with wings	62mm	with	55mm with tube clasp			rotative cap	w/o needle	with rotative cap adaptor	no injection site
62-5F	A62GJ1YLW	latexfree	with air inlet	62mm	with	55mm with tube clasp		with	luer lock	with needle		
62-6F	A62GR2SLO	latexfree	with air inlet	62mm	with	55mm safety	pear-shape		luer lock	w/o needle		safety roller clamp to hold spike isoprene rubber injection site
62-7F	A62GR2YLU	latexfree	with air inlet	62mm	with	55mm safety		with	luer lock	w/o needle	with blue vented cap	safety roller clamp to hold spike
72-1F	A72GJ1LO	latexfree	with air inlet	72mm	with	55mm with tube clasp			luer lock	w/o needle		no injection site
72-2F	A62GJ1XLO	latexfree	with air inlet	62mm	with	55mm with tube clasp			luer lock	w/o needle	anti-reflux valve	no injection site
72-3F	A72GJ1YLU	latexfree	with air inlet	72mm	with	55mm with tube clasp		with	luer lock	w/o needle	with blue vented cap	
72-4F	72GJ1YLO	latexfree	w/o air inlet	72mm	with	55mm with tube clasp		with	luer lock	w/o needle		
1#	50R1NTW	latex	w/o air inlet	50mm		50mm	line-shape		luer tip	with needle		
2# (50-1F)	A50GR1PTW	latex	with air inlet	50mm	with	50mm	pear-shape		luer tip	with needle		
3# (50-3F)	A50GR1PLW	latex	with air inlet	50mm	with	50mm	pear-shape		luer lock	with needle		
4#	A50GJ1YTO	latexfree	with air inlet	50mm	with	55mm with tube clasp		with	luer tip	w/o needle		
5#	A50GJ1YTW	latexfree	with air inlet	50mm	with	55mm with tube clasp		with	luer tip	with needle		
6# (50-2F)	A50GR1YLW	latexfree	with air inlet	50mm	with	50mm		with	luer lock	with needle		
7#	A50GJ1YLO	latexfree	with air inlet	50mm	with	55mm with tube clasp		with	luer lock	w/o needle		

## Configuration of Infusion Set

### Remarks:

1) **Standard Configuration (Economic Configuration):** spike, drip chamber (50mm, 20drops/ml), tubing (length: 150cm, O.D/I.D: 4mm/3mm), regulator (50mm), luer-tip adaptor, with DEHP ingredient

2) **Optional Components to Add:** T-shape spike, spike with wings; air inlet; drip chamber (54mm/62mm/72mm/89mm. 10drops/ml, 60drop/ml); regulator (50mm/55mm with tube clasp/ 55mm with safety hole); solution filtration disk;natural latex /isoprene rubber (latexfree) injection bulb (line-shape/pear-shape); Y-shape injection site(latexfree); needlefree injection site(latexfree); 3-way stopcock injection site(latexfree); extension tubing (O.D/I.D: 4mm/3mm); frosted elastic tubing; male luer lock adaptor; with or without "hypodermic needle, scalp vein set, anti-reflux valve, blue vented cap, micro flow regulator, micro filter" etc; light shielding function etc.; DEHP-free


Model No.	Code of Configuration	Nutral Latex or Latexfree	Spike	Drip Chamber	Solution Filtration Disk	Regulator (Roller Clamp)	Natural Latex Injection Site	Y-shape Injection Site	Adaptor	Needle	Other Components	Remarks
8#	AW54J1YTO	latexfree	with air inlet with wings	54mm	w/o	55mm with tube clasp		with	luer tip	w/o needle		
9#	AW54GR1YLB	latexfree	with air inlet with wings	54mm	with	50mm		with	luer lock	with needle		
10#	J50GR1PTW2	latex	T-shape Spike	50mm	with	50mm	pear-shape		luer tip	with 2 needles		
11#	A50GJ1YTW	latexfree	with air inlet	50mm	with	55mm with tube clasp		with	luer tip	with needle		
12# (54-2F)	AW54GJ1YLO	latexfree	with air inlet with wings	54mm	with	55mm with tube clasp		with	luer lock	w/o needle		
13#	AW54GR1YLO	latexfree	with air inlet with wings	54mm	with	50mm		with	luer lock	w/o needle		
14#	A50GJ1YLB	latexfree	with air inlet	50mm	with	55mm with tube clasp		with	luer lock	with needle		
15#	A50GR1YTO	latexfree	with air inlet	50mm	with	50mm		with	luer tip	w/o needle		
A	A62GR1SLW	latexfree	with air inlet	62mm	with	50mm	pear-shape		luer lock	with needle		isoprene rubber injection bulb
C	A62GJ1SLO	latexfree	with air inlet	62mm	with	55mm with tube clasp	pear-shape		luer lock	w/o needle		isoprene rubber injection bulb
D	A62GR1YLB	latexfree	with air inlet	62mm	with	50mm		with	luer lock	with needle		
E	A72GJ1YLO	latexfree	with air inlet	72mm	with	55mm with tube clasp		with	luer lock	w/o needle		
F (55-1F)	A89GJ1YLO	latexfree	with air inlet	89mm (with ring)	with	55mm with tube clasp		with	luer lock	w/o needle		two-part chamber with rings
G	A62GJ1YLOF	latexfree	with air inlet	62mm	with	55mm with tube clasp		with	luer lock	w/o needle	with micro flow regulator	
S	A50R1NTOV	latex	with air inlet	50mm		50mm	line-shape		luer tip	w/o needle	with scalp vein set	



# Transfusion Set










## Components of Spike, Drip Chamber and Regulator



A- 89mm, for T-shape Spike (w/o air inlet)

B-89mm (w/o air inlet)

C-89mm (with air inlet)

D-89mm(with air inlet, with ring)

E-150mm, double chambers



R1-50mm

J1-55mm with tube clasp

R2-55mm safety



## Configuration of Blood Transfusion Set




**Remarks:**

1) **Standard Configuration (Economic Configuration):** spike, blood drip chamber with filtration net (89mm, 20drops/ml), tubing (length: 150cm, O.D/I.D: 4mm/3mm), regulator (50mm), luer-tip adaptor, with DEHP ingredient


2) **Optional Components to Add:** T-shape spike; air inlet; drip chamber (89mm with ring/150mm); regulator (50mm/55mm with tube clasp/55mm with safety hole); natural latex /isoprene rubber (latexfree) injection bulb (line-shape/pear-shape); Y-shape injection site(latexfree); needlefree injection site(latexfree); extension tubing (O.D/I.D: 4mm/3mm); frosted elastic tubing; male luer lock adaptor; with or without "hypodermic needle, blue vented cap" etc; light shielding function etc.; DEHP-free

Model No.	Code of Configuration	Nutral Latex or Latexfree	Spike	Drip Chamber	Regulator (Roller Clamp)	Latex Injection Site	Y-shape Injection Site	Adaptor	Needle	Other Components	Remarks
B1F (T2)	A89J1PLW	latex	with air inlet	89mm	55mm with tube clasp	pear-shape		luer lock	with needle		
B2F	A89HJ1YLU	latexfree	with air inlet	89mm with ring	55mm with tube clasp		with	luer lock	w/o needle	blue vented cap	
B3F	AJ89R1LO	latexfree	with air inlet	89mm	50mm			luer lock	w/o needle	2 spikes, 2 regulartors	no injection site
B4F	89J1YLU	latexfree	w/o air inlet	89mm	55mm with tube clasp		with	luer lock	with needle		
B5F	150J1YLU	latexfree	w/o air inlet	150mm	55mm with tube clasp		with	luer lock	with needle		double-chamber drip chamber, total 150mm
B6F	A150J1LOC	latexfree	with air inlet	150mm	55mm with tube clasp			luer lock	w/o needle	needlefree injection site	double-chamber drip chamber, total 150mm
T1	A89J1PTW	latex	with air inlet	89mm	55mm with tube clasp	pear-shape		luer tip	with needle		
T2 (B1F)	A89J1PLW	latex	with air inlet	89mm	55mm with tube clasp	pear-shape		luer lock	with needle		
T3	A89HJ1PLO	latex	with air inlet	89mm with ring	55mm with tube clasp	pear-shape		luer lock	w/o needle		drip chamber with ring
T4	A89HJ1YLO	latexfree	with air inlet	89mm with ring	55mm with tube clasp		with	luer lock	w/o needle		drip chamber with ring
T5	A89R1YLO	latexfree	with air inlet	89mm	50mm		with	luer lock	w/o needle		
T6	A150J1YLU	latexfree	with air inlet	150mm	55mm with tube clasp		with	luer lock	with needle		double-chamber drip chamber, total 150mm
T7	J89R1PTO	latex	T-shape w/o air inlet	89mm	50mm	pear-shape		luer tip	w/o needle	2 spikes, 3 regulartors	

# Extension Tube



## I.V Cannula, Dental Needle & Blood Collection Needle for Vacuum Tube Use



I.V Cannula  
model A1: pen model, shell package  
model A2: pen model, blister package  
model B: with wing  
model C: with injection valve



Dental Needle



Blood Collection Needle for Vacuum Tube Use



## Disposable Scalp Vein Set


- Connecting End: Luer Tip / Luer Tip with cap / Luer Lock with screw cap; • Wings: **A**-Single Wing / **B**-Butterfly Wings;
- Individual Package: PE bag / soft blister bag / hard blister pack; • Packing: 100pcs/Papercard box, 2000pcs/ctn; • Sterile: EO Sterile
- Specifications:

Size		Color	Regular Needle Length					Regular Length
OD (mm)	Gauge		1/2"	5/8"	3/4"	1"	1 1/8"	
0.40	27G	Grey	•		•			300mm (12")
0.45	26G	Brown		•	•			
0.50	25G	Orange			•			
0.55	24G	Purple			•			
0.60	23G	Blue			•	•	•	
0.70	22G	Black			•	•	•	
0.80	21G	Green			•	•	•	
0.90	20G	Yellow			•		•	
1.00	19G	Beige			•		•	
1.20	18G	Pink			•		•	



## Disposable Arterial Venous Fistula Needle (A.V. Needle)

- with / without back eye
- Fixed / Rotated Wings
- Sizes: 14G, 15G, 16G, 17G
- Tube length: 30cm
- Individual Package: PE bag  
Paper-PE Peel Pack
- Packing: 50pcs/box, 1000pcs/ctn
- Sterile: EO Sterile



## Disposable Dental Needle

- **Model A:** 25G (Orange): Length (mm): 13, 21, 25, 32, 38, 42; 27G (Brown): Length (mm): 13, 21, 25, 32, 38, 42; 30G (Yellow): Length (mm): 13, 17, 21, 25, 32, 38
- **Model B:** 25G/24G/23G/22G/21G/20G/18G
- Individual Package: Sticker Seal • Packing: 100pcs/box, 5000pcs/ctn • Sterile: EO Sterile

